

The module for automating your dosing tasks



HandyStep[®] electronic

with RS-232 interface

F I R S T C L A S S · B R A N D

Specially designed for integration into automated systems, the HandyStep[®] electronic with interface forms the main component for solving complex dispensing tasks.

- **External controllable motor dispenser** as dispensing unit for individual automation
- **Dispensing, sequential dispensing and pipetting** easily in remote mode
- **Aspiration and dispensing speed** can be adjusted independently
- **Variable volume adjustment** from 1.0 µl to 50 ml
- **Use of high-quality BRAND[®] Precision Dispenser Tips** in 10 different sizes and available sterile if necessary
- **CE-IVD-compliant**



The BRAND HandyStep® electronic hand dispenser with RS-232 interface is best suited for streamlining time-consuming dosing and filling processes and is also known for its high reliability. The quality of the servo motor unit guarantees reproducible dosing results.

Precision Dispenser Tips

The use of BRAND PD-Tips as disposable dosing units (cylinder: PP / Piston: PE-HD, 0.1



ml: LCP) allows a simple and quick change in the production batch or of the dosing medium. Only the original BRAND PD-Tips are ideally suited for remotely operating the HandyStep® electronic. These operate on the positive displacement principle and are suitable for dispensing highly viscous or high vapor pressure fluids and foaming media. The error limits of ISO 8655-5 are complied with.

ml: LCP) allows a simple and quick change in the production batch or of the dosing medium.

Only the original BRAND PD-Tips are ideally suited for remotely operating the HandyStep® electronic. These operate on the positive displacement principle and are suitable for dispensing highly viscous or high vapor pressure fluids and foaming media. The error limits of ISO 8655-5 are complied with.

Remote mode functions:

The following customized dispensing parameters can be programmed and controlled from the primary host (host computer):

- Volume setting from 1.0 µl to 50 ml
- Also, dispensing can be programmed for various volumes using only a single PD-Tip filling (sequential dosing), e.g. first step 10 µl, second step 4 ml, etc.
- Speed for aspiration and dispensing can be adjusted independently. Three speed levels are available.
- Two modes available: Multi-volume dispensing or single-volume pipetting

Ordering Information



HandyStep® electronic with RS-232 interface

with socket, incl. NiMH battery pack, charger and power supply

Cat. No. 7050 80



Also comes with the following:
Software for WINDOWS® 95 and later operating systems, Library HSeBRAND.DLL, **2 m interface cable**

Software with source files for transferring to other operating systems
- Delivery on request -

Detailed information about the HandyStep® electronic (standard unit without interface) and PD-Tips on the Internet at www.brand.de.

HandyStep® and BRAND® are trademarks of BRAND GMBH + CO KG, Germany.

Other reproduced brands are the property of the respective owner.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG · P.O. Box 11 55 · 97861 Wertheim · Germany

Tel.: +49 9342 808-0 · Fax: +49 9342 808-98000 · E-Mail: info@brand.de · Internet: www.brand.de

